



Avazzia Star Product Description

This product is a hand-held, battery-operated, microcurrent electro-stimulation medical device. The Avazzia Star is designed for at-home use and offers easy-to-use features, including four frequencies that allow patients to administer micro-current pain relief at home.

Research has shown microcurrent offers long-term pain relief without drugs or surgery. It is not habit-forming and non-addictive.



Avazzia Star

PRODUCT SPECIFICATIONS	
SIZE AND WEIGHT	Size: 117mm X 64mm X 33mm. Weight (device only): Approx. 103 grams (3.6 oz). Weight (with batteries): Approx. 150 grams (5.3 oz)
FREQUENCY, HZ	28 Hz to 222 Hz
POWER SOURCE	2 AA 1.5V DC batteries
OUTPUT	Output Voltage: 20 V to 500 V Output Current: 2 mA to 90 mA Current Density: 0.1 to 966 uA/cm ² Power Density: 0.1 to 1.2 mW/cm ²
PULSES	Pulse Duration: less than 2 ms
WAVEFORM	Pulsed, damped, biphasic sinusoidal
AUTOMATIC SHUT-OFF	60 minutes

PRODUCT INFORMATION	
MODES OF OPERATION	AS 77-3, AS 77-8, AS 148-2, AS STIM
VISUAL INDICATIONS	Mode selection and power intensity indicating LEDs
PRODUCTION	Manufactured by Avazzia in U.S.
AVAILABILITY	Medical device available for general purchase. Visit www.avazzia.com for information.
ACCESSORIES (SOLD SEPARATELY)	Y-electrode, adhesive conductive pads, brush electrode, pencil electrode
CLEARANCES	ISO 13485 international certification CE approved, CB certified HealthCanada certified
WARNINGS	Not to be used by those with pacemaker or other implantable electrical device. Not to be used by pregnant or nursing women. Do not use directly over carotid artery or cause current to run transcranially. Do not use directly on open wounds.